



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PRKAB1 Antibody [#C12821]

Catalog Number: C12821

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y478

Other Names: 5'-AMP-activated protein kinase subunit beta-1, AMPK beta-1 chain, AMPKb, PRKAB1, AMPK

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human PRKAB1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PRKAB1 antibody detects endogenous levels of total PRKAB1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Stapleton D., FEBS Lett. 409:452-456(1997).

Scherer S.E., Nature 440:346-351(2006).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

