



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## BARD1 Antibody [#C13026]

**Catalog Number:** C13026

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q99728

**Other Names:** BRCA1-associated RING domain protein 1, BARD-1, BARD1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human BARD1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** BARD1 antibody detects endogenous levels of total BARD1 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Wu L.C., Nat. Genet. 14:430-440(1996).

Thai T.H., Hum. Mol. Genet. 7:195-202(1998).

Hillier L.W., Nature 434:724-731(2005).

