



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## DCC Antibody [#C13040]

**Catalog Number:** C13040

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P43146

**Other Names:** colorectal cancer suppressor, deleted in colorectal carcinoma, netrin receptor DCC, tumor suppressor protein DCC

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human DCC.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** DCC antibody detects endogenous levels of total DCC protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 100%, Positives = 100%);  
Rat (Identities = 100%, Positives = 100%)

**Applications:** IF: 1:100~1:500 ELISA: 1:20000

**References:**

Hedrick L., Genes Dev. 8:1174-1183(1994).

Cho K.R., Genomics 19:525-531(1994).

Sjoeblom T., Science 314:268-274(2006).

