



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## ELOA1 Antibody [#C13047]

**Catalog Number:** C13047

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q14241

**Other Names:** Transcription elongation factor B polypeptide 3, RNA polymerase II transcription factor SIII subunit A1, SIII p110, Elongin-A, EloA, Elongin 110 kDa subunit, TCEB3, MSTP059

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

ELOA1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ELOA1 antibody detects endogenous levels of total ELOA1 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:1000

**References:**

Aso T., Gene 168:277-278(1996).

Gregory S.G., Nature 441:315-321(2006).


Olsen J.V., Cell 127:635-648(2006).



Peptide -

+

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using ELOA1 antibody.



Western blot analysis of extracts from Jurkat cells and HUVEC cells, using ELOA1 antibody.