



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ELOA3 Antibody [#C13049]

Catalog Number: C13049

Concentration: 1mg/ml

Swiss-Prot No.: Q8NG57

Other Names: RNA polymerase II transcription factor SIII subunit A3, Elongin-A3, EloA3, Transcription elongation factor B polypeptide 3C, TCEB3C, TCEB3L2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human ELOA3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ELOA3 antibody detects endogenous levels of total ELOA3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Yamazaki K., J. Biol. Chem. 277:26444-26451(2002).

