



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## CDK5RAP3 Antibody [#C12192]

**Catalog Number:** C12192

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q96JB5

**Other Names:** CDK5 regulatory subunit-associated protein 3, CDK5 activator-binding protein C53, Protein HSF-27, IC53, MSTP016, OK/SW-cl.114, PP1553, CK5P3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

CDK5RAP3.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CDK5RAP3 antibody detects endogenous levels of total CDK5RAP3 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Chen J., Biochem. Biophys. Res. Commun. 294:161-166(2002).

Xie Y.H., Cell Res. 13:83-91(2003).

Ota T., Nat. Genet. 36:40-45(2004).

