



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CMKLR1 Antibody [#C12193]

Catalog Number: C12193

Concentration: 1mg/ml

Swiss-Prot No.: Q99788

Other Names: Chemokine receptor-like 1, G-protein coupled receptor DEZ, G-protein coupled receptor ChemR23, CHEMR23, DEZ, CML1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CMKLR1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CMKLR1 antibody detects endogenous levels of total CMKLR1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Methner A., Biochem. Biophys. Res. Commun. 233:336-342(1997).

Samson M., Eur. J. Immunol. 28:1689-1700(1998).

Hillman R.T., Genome Biol. 5:RESEARCH008.1-RESEARCH008.16(2004).

