



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Collagen IX α 3 Antibody [#C12212]

Catalog Number: C12212

Concentration: 1mg/ml

Swiss-Prot No.: Q14050

Other Names: Collagen alpha-3(IX) chain , CO9A3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Collagen IX α 3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Collagen IX α 3 antibody detects endogenous levels of total Collagen IX α 3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:5000

References:

Brewton R.G., Genomics 30:329-336(1995).


Paasilta P., J. Biol. Chem. 274:22469-22475(1999).

Deloukas P., Nature 414:865-871(2001).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human tonsil tissue using Collagen IX α 3 antibody.



Peptide - +

Immunofluorescence analysis of NIH/3T3 cells, using Collagen IX α 3 antibody.