



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## LAMB2 Antibody [#C13070]

**Catalog Number:** C13070

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P55268

**Other Names:** laminin B1s; laminin beta-2; laminin S;  
LAMS; S-laminin

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human LAMB2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** LAMB2 antibody detects endogenous levels of total LAMB2 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Wewer U.M., Genomics 24:243-252(1994).

Iivanainen A., Matrix Biol. 14:489-497(1995).

Zenker M., Hum. Mol. Genet. 13:2625-2632(2004).

