



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LAMC1 Antibody [#C13072]

Catalog Number: C13072

Concentration: 1mg/ml

Swiss-Prot No.: P11047

Other Names: Laminin subunit gamma-1 [Precursor],
Laminin B2 chain, LAMC1, LAMB2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human LAMC1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LAMC1 antibody detects endogenous levels of total LAMC1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Kallunki T., J. Biol. Chem. 266:221-228(1991).

Pikkarainen T., J. Biol. Chem. 263:6751-6758(1988).

Fukushima Y., Cytogenet. Cell Genet. 48:137-141(1988).

