



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## MLH3 Antibody [#C13087]

**Catalog Number:** C13087

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9UHC1

**Other Names:** DNA mismatch repair protein Mlh3, MutL protein homolog 3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human MLH3.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** MLH3 antibody detects endogenous levels of total MLH3 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Lipkin S.M., Nat. Genet. 24:27-35(2000).

Kondo E., Nucleic Acids Res. 29:1695-1702(2001).

Lipkin S.M., Hum. Mutat. 17:389-396(2001).

