



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MSH6 Antibody [#C13091]

Catalog Number: C13091

Concentration: 1mg/ml

Swiss-Prot No.: P52701

Other Names: DNA mismatch repair protein Msh6, MutS-alpha 160 kDa subunit, G/T mismatch-binding protein, GTBP, GTMBP, p160, MSH6, GTBP

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MSH6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MSH6 antibody detects endogenous levels of total MSH6 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Acharya S., Proc. Natl. Acad. Sci. U.S.A. 93:13629-13634(1996).

Shiwaku H.O., DNA Res. 4:359-362(1997).

Palombo F., Science 268:1912-1914(1995).

