



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## NBN Antibody [#C13093]

**Catalog Number:** C13093

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O60934

**Other Names:** Nibrin, Nijmegen breakage syndrome protein 1, Cell cycle regulatory protein p95, NBN, NBS, NBS1, P95

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human NBN.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** NBN antibody detects endogenous levels of total NBN protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Varon R., Cell 93:467-476(1998).

Carney J.P., Cell 93:477-486(1998).

Matsuura S., Nat. Genet. 19:179-181(1998).

