



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PARP2 Antibody [#C13097]

Catalog Number: C13097

Concentration: 1mg/ml

Swiss-Prot No.: Q9UGN5

Other Names: Poly [ADP-ribose] polymerase 2; PARP-2; NAD(+) ADP-ribosyltransferase 2; Poly[ADP-ribose] synthetase 2; pADPRT-2; hPARP-2; ADPRT2; ADPRTL2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human PARP2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PARP2 antibody detects endogenous levels of total PARP2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)
Mouse (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:40000

References:

Sugano S., Nat. Genet. 36:40-45(2004).

de Murcia G.M., J. Biol. Chem. 277:23028-23036(2002).

Bennett K.L., Submitted (OCT-2008).

