



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PARP3 Antibody [#C13098]

Catalog Number: C13098

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y6F1

Other Names: ADP-ribosyltransferase 3, ADPRT3, ADPRTL2, ADPRTL3, EC 2.4.2.30, IRT1, NAD+ ADP-ribosyltransferase 3, PARP-3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human PARP3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PARP3 antibody detects endogenous levels of total PARP3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:1000

References:

Johansson M., Genomics 57:442-445(1999).

Augustin A., J. Cell Sci. 116:1551-1562(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue using PARP3 antibody.