



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RAD21 Antibody [#C13110]

Catalog Number: C13110

Concentration: 1mg/ml

Swiss-Prot No.: O60216

Other Names: Double-strand-break repair protein rad21 homolog, hHR21, Nuclear matrix protein 1, NXP-1, SCC1 homolog, RAD21, HR21, KIAA0078, NXP1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RAD21.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RAD21 antibody detects endogenous levels of total RAD21 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

McKay M.J., Genomics 36:305-315(1996).

Sadano H., Biochem. Biophys. Res. Commun. 267:418-422(2000).

Nomura N., DNA Res. 1:223-229(1994).

