



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

FAF1 Antibody [#C13054]

Catalog Number: C13054

Concentration: 1mg/ml

Swiss-Prot No.: Q9UNN5

Other Names: CGI-03; FAS-associated factor 1; Hfaf1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human FAF1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FAF1 antibody detects endogenous levels of total FAF1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Ryu S.-W., Biochem. Biophys. Res. Commun. 262:388-394(1999).

Lai C.-H., Genome Res. 10:703-713(2000).

Bechtel S., BMC Genomics 8:399-399(2007).

