



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## LAMA1 Antibody [#C13064]

**Catalog Number:** C13064

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P25391

**Other Names:** LAMA, laminin A, laminin alpha-1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human LAMA1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** LAMA1 antibody detects endogenous levels of total LAMA1 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%)

**Applications:** IF: 1:100~1:500 ELISA: 1:40000

**References:**

Haaparanta T., Matrix 11:151-160(1991).

Nissinen M., Biochem. J. 276:369-379(1991).

Olsen D., Lab. Invest. 60:772-782(1989).

