



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RAD50 Antibody [#C13111]

Catalog Number: C13111

Concentration: 1mg/ml

Swiss-Prot No.: Q92878

Other Names: DNA repair protein RAD50, EC 3.6.-.-, hRAD50

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RAD50.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RAD50 antibody detects endogenous levels of total RAD50 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Dolganov G.M., Mol. Cell. Biol. 16:4832-4841(1996).

Kim K.K., Gene 235:59-67(1999).

Carney J.P., Cell 93:477-486(1998).

