



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Serpin B5 Antibody [#C13117]

Catalog Number: C13117

Concentration: 1mg/ml

Swiss-Prot No.: P36952

Other Names: Serpin B5 [Precursor], Maspin, Protease inhibitor 5, SERPINB5, PI5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Serpin B5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Serpin B5 antibody detects endogenous levels of total Serpin B5 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Zou Z., Science 263:526-529(1994).

Pemberton P.A., J. Biol. Chem. 270:15832-15837(1995).

Law R.H., J. Biol. Chem. 280:22356-22364(2005).

