



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TIPI1 Antibody [#C13126]

Catalog Number: C13126

Concentration: 1mg/ml

Swiss-Prot No.: Q9HCN2

Other Names: p53-regulated apoptosis-inducing protein 1, P53aip1, TP53AIP1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TIPI1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TIPI1 antibody detects endogenous levels of total TIPI1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:1000

References:

Oda K., Cell 102:849-862(2000).


Zhou G., Submitted (MAY-2003) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human thymus gland tissue using TIPI1 antibody.



Peptide - +

Immunofluorescence analysis of A549 cells, using TIPI1 antibody.