



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

POU4F3 Antibody [#C13491]

Catalog Number: C13491

Concentration: 1mg/ml

Swiss-Prot No.: Q15319

Other Names: POU domain class 4 transcription factor 3;
Brain-specific homeobox/POU domain protein 3C; Brn-3C;
POU4F3; BRN3C

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human POU4F3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: POU4F3 antibody detects endogenous levels of total POU4F3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Xiang M., J. Neurosci. 15:4762-4785(1995).

Vahava O., Science 279:1950-1954(1998).

Collin R.W.J., Hum. Mutat. 29:545-554(2008).

