



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AGPAT5 Antibody [#C14006]

Catalog Number: C14006

Concentration: 1mg/ml

Swiss-Prot No.: Q9NUQ2

Other Names: 1-acyl-sn-glycerol-3-phosphate acyltransferase epsilon; 1-AGP acyltransferase 5; 1-AGPAT 5; 1-AGPAT5; EC 2.3.1.51; PLCE

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AGPAT5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AGPAT5 antibody detects endogenous levels of total AGPAT5 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Wiemann S., Genome Res. 11:422-435(2001).

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).



For Research Use Only