



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AGPAT3 Antibody [#C14008]

Catalog Number: C14008

Concentration: 1mg/ml

Swiss-Prot No.: Q9NRZ7

Other Names: 1-acylglycerol-3-phosphate
o-acyltransferase 3; agpat3; agpat3 protein;
lpaat-gamma1; lysophosphatidic acid
acyltransferase-gamma1; lysophosphatidic acid
acyltransferase-gamma2; mgc4604; plcc

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human

AGPAT3.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: AGPAT3 antibody detects endogenous levels
of total AGPAT3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 96%, Positives = 100%);
Rat (Identities = 96%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Leung D.W.; Front. Biosci. 6:D944-D953(2001).

Clark H.F., Genome Res. 13:2265-2270(2003).

Ota T., Nat. Genet. 36:40-45(2004).

