



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## PSMC6 Antibody [#C14017]

**Catalog Number:** C14017

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P62333

**Other Names:** 26S protease regulatory subunit S10B; proteasome 26S subunit ATPase 6; proteasome subunit p42; PRS10; SUG2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human PSMC6.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** PSMC6 antibody detects endogenous levels of total PSMC6 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%)  
Mouse (Identities = 100%, Positives = 100%)

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Fujiwara T., FEBS Lett. 387:184-188(1996).

Kalnine N., Submitted (MAY-2003) to the EMBL/GenBank/DDBJ databases.

Demartino G.N., J. Biol. Chem. 271:3112-3118(1996).

