



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PSMD11 Antibody [#C14020]

Catalog Number: C14020

Concentration: 1mg/ml

Swiss-Prot No.: O00231

Other Names: 26S proteasome non-ATPase regulatory subunit 11; 26S proteasome regulatory subunit p44.5; 26S proteasome regulatory subunit S9; PSD11

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human PSMD11.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PSMD11 antibody detects endogenous levels of total PSMD11 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Hoffman L., FEBS Lett. 404:179-184(1997).

Saito A., Gene 203:241-250(1997).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

