



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

MRPS12 Antibody [#C14033]

Catalog Number: C14033

Concentration: 1mg/ml

Swiss-Prot No.: O15235

Other Names: 28S ribosomal protein S12 mitochondrial; S12mt; MRP-S12; MT-RPS12

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MRPS12.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MRPS12 antibody detects endogenous levels of total MRPS12 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 86%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Shah Z.H., Gene 204:55-62(1997).

Johnson D.F., Genomics 52:363-368(1998).

Ebert L., Submitted (MAY-2004) to the EMBL/GenBank/DDBJ databases.

