



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MRPS9 Antibody [#C14045]

Catalog Number: C14045

Concentration: 1mg/ml

Swiss-Prot No.: P82933

Other Names: 28S ribosomal protein S9, mitochondrial;
MRP-S9; RPMS9; RT09; S9mt

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MRPS9.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MRPS9 antibody detects endogenous levels of total MRPS9 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)
Mouse (Identities = 98%, Positives = 194%);
Rat (Identities = 94%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

The MGC Project Team; Genome Res. 14:2121-2127(2004).
Koc E.C., J. Biol. Chem. 276:19363-19374(2001).

