



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MRPL10 Antibody [#C14054]

Catalog Number: C14054

Concentration: 1mg/ml

Swiss-Prot No.: Q7Z7H8

Other Names: l10mt; mgc17973; mitochondrial ribosomal protein l10; mrp-l10; mrp-l8; mrpl10; mrpl8; rm10; rpml8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MRPL10.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MRPL10 antibody detects endogenous levels of total MRPL10 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Kenmochi N., Genomics 77:65-70(2001).

Koc E.C., J. Biol. Chem. 276:43958-43969(2001).

