



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MRPL12 Antibody [#C14056]

Catalog Number: C14056

Concentration: 1mg/ml

Swiss-Prot No.: P52815

Other Names: 39S ribosomal protein L12, mitochondrial; 5c5-2; L12mt; MRP-L12; RM12; RPML12

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MRPL12.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MRPL12 antibody detects endogenous levels of total MRPL12 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)
Mouse (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Marty L., J. Biol. Chem. 271:11468-11476(1996).

Liu J., Submitted (APR-1998) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

