



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## MRPL14 Antibody [#C14058]

**Catalog Number:** C14058

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q6P1L8

**Other Names:** 39S ribosomal protein L14, mitochondrial; L14mt; MRP-L14; MRP-L32; RM14; RPML32

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human MRPL14.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** MRPL14 antibody detects endogenous levels of total MRPL14 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Mungall A.J., Nature 425:805-811(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Kenmochi N., Genomics 77:65-70(2001).

