



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RPS8 Antibody [#C14126]

Catalog Number: C14126

Concentration: 1mg/ml

Swiss-Prot No.: P62241

Other Names: 40S ribosomal protein S8; OK/SW-cl.83; ribosomal protein S8; RS8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RPS8.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RPS8 antibody detects endogenous levels of total RPS8 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Davies B., Genomics 15:68-75(1993).

Shichijo S., Submitted (MAY-2001) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

