



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RPS9 Antibody [#C14127]

Catalog Number: C14127

Concentration: 1mg/ml

Swiss-Prot No.: P46781

Other Names: 40S ribosomal protein S9; ribosomal protein S9; RS9

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RPS9.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RPS9 antibody detects endogenous levels of total RPS9 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Frigerio J.-M., *Biochim. Biophys. Acta* 1262:64-68(1995).

Yoshihama M., *Genome Res.* 12:379-390(2002).

Vladimirov S.N., *Eur. J. Biochem.* 239:144-149(1996).



Western blot analysis of extracts from HepG2 cells, COLO cells, HT-29 cells, HeLa cells and HUVEC cells, using RPS9 antibody.



Immunofluorescence analysis of A549 cells, using RPS9 antibody.