



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RPL14 Antibody [#C14151]

Catalog Number: C14151

Concentration: 1mg/ml

Swiss-Prot No.: P50914

Other Names: CAG-ISL-7; CTG-B33; hRL14; L14; MGC88594; ribosomal protein L14; RL14

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RPL14.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RPL14 antibody detects endogenous levels of total RPL14 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 92%, Positives = 100%);
Rat (Identities = 92%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Aoki M., Diabetes 45:157-164(1996).

Tanaka M., Biochem. Biophys. Res. Commun. 243:531-537(1998).

Yoshihama M., Genome Res. 12:379-390(2002).

