



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RPL26L Antibody [#C14163]

Catalog Number: C14163

Concentration: 1mg/ml

Swiss-Prot No.: Q9UNX3

Other Names: 60S ribosomal protein L26-like 1;
RPL26L1; RPL26P1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human RPL26L.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RPL26L antibody detects endogenous levels of total RPL26L protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Zhang Q.-H., Genome Res. 10:1546-1560(2000).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

