



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RPL28 Antibody [#C14166]

Catalog Number: C14166

Concentration: 1mg/ml

Swiss-Prot No.: P46779

Other Names: FLJ43307; ribosomal protein L28; RL28

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RPL28.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RPL28 antibody detects endogenous levels of total RPL28 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Frigerio J.-M., *Biochim. Biophys. Acta* 1262:64-68(1995).
Ebert L., Submitted (JUN-2004) to the EMBL/GenBank/DDBJ databases.
The MGC Project Team; *Genome Res.* 14:2121-2127(2004).

