



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RPL35 Antibody [#C14173]

Catalog Number: C14173

Concentration: 1mg/ml

Swiss-Prot No.: P42766

Other Names: 60S ribosomal protein L35; ribosomal protein L35; RL35

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RPL35.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RPL35 antibody detects endogenous levels of total RPL35 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Patel S.K., Submitted (JUL-1994) to the EMBL/GenBank/DDBJ databases

Ebert L., Submitted (JUN-2004) to the EMBL/GenBank/DDBJ databases.

Humphray S.J., Nature 429:369-374(2004).

