



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## RPL7 Antibody [#C14189]

**Catalog Number:** C14189

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P18124

**Other Names:** 60S ribosomal protein L7; RL7; RPL7

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human RPL7.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** RPL7 antibody detects endogenous levels of total RPL7 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Hemmerich P., Nucleic Acids Res. 21:223-231(1993).

Seshadri T., J. Biol. Chem. 268:18474-18480(1993).

Schneider R., Submitted (MAY-1990) to the EMBL/GenBank/DDBJ databases.

