



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ABHD14A Antibody [#C14212]

Catalog Number: C14212

Concentration: 1mg/ml

Swiss-Prot No.: Q9BUJ0

Other Names: Abhydrolase domain-containing protein 14A; EC 3.-.-.-

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ABHD14A.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ABHD14A antibody detects endogenous levels of total ABHD14A protein.

Reactivity: Human (Identities = 100%, Positives = 100%)
Mouse (Identities = 92%, Positives = 100%);
Rat (Identities = 92%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Clark H.F., Genome Res. 13:2265-2270(2003).

The German cDNA consortium; Submitted (SEP-2004) to the EMBL/GenBank/DDBJ databases.

Genome Res. 14:2121-2127(2004).

