



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ACTR3 Antibody [#C14244]

Catalog Number: C14244

Concentration: 1mg/ml

Swiss-Prot No.: P61158

Other Names: Actin-2; Actin-like protein 3; Actin-related protein 3; ACTR3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ACTR3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ACTR3 antibody detects endogenous levels of total ACTR3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Welch M.D., J. Cell Biol. 138:375-384(1997).

Welch M.D., Nature 385:265-269(1997).

Gournier H., Mol. Cell 8:1041-1052(2001).

