



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## ACAD10 Antibody [#C14253]

**Catalog Number:** C14253

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q6JQN1

**Other Names:** Acyl-CoA dehydrogenase family member;  
10ACAD-10

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ACAD10.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ACAD10 antibody detects endogenous levels of total ACAD10 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%)


**Applications:** IHC: 1:50~1:100 ELISA: 1:40000

**References:**

Ye X., Mol. Biol. Rep. 31:191-195(2004).

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human heart tissue using ACAD10 antibody.