



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ACOT2 Antibody [#C14271]

Catalog Number: C14271

Concentration: 1mg/ml

Swiss-Prot No.: P49753

Other Names: acyl-coenzyme A thioesterase 2; CTE-Ia;
EC 3.1.2.2; peroxisomal acyl-coenzyme A thioester
hydrolase 2a; peroxisomal long-chain acyl-coA
thioesterase 2; ZAP128

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human

ACOT2.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: ACOT2 antibody detects endogenous levels
of total ACOT2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Sherrington R., Nature 375:754-760(1995).

Hunt M.C., FASEB J. 20:1855-1864(2006).

Ota T., Nat. Genet. 36:40-45(2004).

