



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADAMTS18 Antibody [#C14310]

Catalog Number: C14310

Concentration: 1mg/ml

Swiss-Prot No.: Q8TE60

Other Names: A disintegrin and metalloproteinase with thrombospondin motifs 18, ADAM-TS 18, ADAM-TS18, ADAMTS-18EC=3.4.24.-

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ADAMTS18.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADAMTS18 antibody detects endogenous levels of total ADAMTS18 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Cal S., Gene 283:49-62(2002) [PubMed: 11867212].

Martin J., Nature 432:988-994(2004) [PubMed: 15616553].

Ota T., Nat. Genet. 36:40-45(2004) [PubMed: 14702039].



For Research Use Only