



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SLC25A6 Antibody [#C14334]

Catalog Number: C14334

Concentration: 1mg/ml

Swiss-Prot No.: P12236

Other Names: adenine nucleotide translocator 2; adenine nucleotide translocator 3; ADP,ATP carrier 3; ADP,ATP carrier, isoform T2; ADP,ATP carrier, liver isoform T2; ADP/ATP translocase 3; ADP/ATP translocator of liver; ADT3; ANT 3; ANT3; ANT3Y; solute carrier family 25 member 6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

SLC25A6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SLC25A6 antibody detects endogenous levels of total SLC25A6 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Cozens A.L., J. Mol. Biol. 206:261-280(1989).

Zhou J., Submitted (JUL-2000) to the EMBL/GenBank/DDBJ databases.

Houldsworth J., Proc. Natl. Acad. Sci. U.S.A. 85:377-381(1988).

