



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

BST1 Antibody [#C14339]

Catalog Number: C14339

Concentration: 1mg/ml

Swiss-Prot No.: Q10588

Other Names: ADP-ribosyl cyclase 2; bone marrow stromal antigen 1; BST-1; cADPr hydrolase 2; CD157 antigen; cyclic ADP-ribose hydrolase 2; EC 3.2.2.5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human BST1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: BST1 antibody detects endogenous levels of total BST1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Kaisho T., Proc. Natl. Acad. Sci. U.S.A. 91:5325-5329(1994).

Kalnine N., Submitted (OCT-2004) to the EMBL/GenBank/DDBJ databases.

Hillier L.W., Nature 434:724-731(2005).

