



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARF6 Antibody [#C14344]

Catalog Number: C14344

Concentration: 1mg/ml

Swiss-Prot No.: P62330

Other Names: ADP-ribosylation factor 6; ARF6;

DKFZp564M0264

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human ARF6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: ARF6 antibody detects endogenous levels of total ARF6 protein.

Reactivity: Human; Mouse; Rat


Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Tsuchiya M., J. Biol. Chem. 266:2772-2777(1991).

Mao M., Proc. Natl. Acad. Sci. U.S.A. 95:8175-8180(1998).

Puhl H.L. III, Submitted (MAR-2002) to the EMBL/GenBank/DDBJ databases.



For Research Use Only