



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARFGAP3 Antibody [#C14346]

Catalog Number: C14346

Concentration: 1mg/ml

Swiss-Prot No.: Q9NP61

Other Names: ADP-ribosylation factor GTPase activating protein 1; ADP-ribosylation factor GTPase activating protein 3; ADP-ribosylation factor GTPase-activating protein 3; ARF GAP 3; ARFG3; ARFGAP1; ARFGAP3; FLJ45618

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human

ARFGAP3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARFGAP3 antibody detects endogenous levels of total ARFGAP3 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Zhang C., Genomics 63:400-408(2000).

Collins J.E., Genome Res. 13:27-36(2003).

Collins J.E., Genome Biol. 5:R84.1-R84.11(2004).



For Research Use Only