



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LPCAT2 Antibody [#C14281]

Catalog Number: C14281

Concentration: 1mg/ml

Swiss-Prot No.: Q7L5N7

Other Names: acyltransferase-like 1; EC 2.3.1.-; PCAT2;
PCATA

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human LPCAT2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: LPCAT2 antibody detects endogenous levels of total LPCAT2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)
Mouse (Identities = 94%, Positives = 100%);
Rat (Identities = 94%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

The German cDNA consortium; Submitted (AUG-2003) to the EMBL/GenBank/DDBJ databases.

