



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADH7 Antibody [#C14383]

Catalog Number: C14383

Concentration: 1mg/ml

Swiss-Prot No.: P40394

Other Names: ADH-4; alcohol dehydrogenase class 4 mu/sigma; alcohol dehydrogenase class IV mu/sigma; EC 1.1.1.1; gastric alcohol dehydrogenase; retinol dehydrogenase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ADH7.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADH7 antibody detects endogenous levels of total ADH7 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Satre M.A., J. Biol. Chem. 269:15606-15612(1994).

Zgombic-Knight M., J. Biol. Chem. 270:4305-4311(1995).

Farres J., Eur. J. Biochem. 224:549-557(1994).

