



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

ALDH1B1 Antibody [#C14388]

Catalog Number: C14388

Concentration: 1mg/ml

Swiss-Prot No.: P30837

Other Names: acetaldehyde dehydrogenase 5; AL1B1; aldehyde dehydrogenase X, mitochondrial; ALDH class 2; ALDH5; ALDHX; EC 1.2.1.3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ALDH1B1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ALDH1B1 antibody detects endogenous levels of total ALDH1B1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Hsu L.C., J. Biol. Chem. 266:12257-12265(1991).

Kalnine N., Submitted (MAY-2003) to the EMBL/GenBank/DDBJ databases.

Humphray S.J., Nature 429:369-374(2004).

